


# Onx Low Anticoagulation Study Results



**2001** US FDA approves On-X Aortic Valve

**2002** US FDA approves On-X Mitral Valve

US FDA approves IDE lowered anticoagulation clinical for the On-X heart valve.

**2006** Prospective Randomized On-X Valve Anticoagulation Clinical Trial (PROACT) is initiated.

**2009** Enrollment complete in the high-risk aortic patient group of the Prospective Randomized On-X Valve Anticoagulation Clinical Trial (PROACT).

**2011** US FDA approves On-X Ascending Aortic Prosthesis

**2013** US FDA approves Chord-X™ ePTFE suture for mitral valve repair and replacement

**2013** US FDA approves On-X aortic heart valve with anatomic sewing cuff

**2013** Enrollment complete in the low-risk aortic patient group and the mitral patient group of the Prospective Randomized On-X Valve Anticoagulation Clinical Trial (PROACT)

**US FDA approves On-X aortic heart valves with INR 1.5–2.0 after 3 months standard**

**2015** anticoagulation, based on results from the high risk AVR group of the Prospective Randomized On-X Valve Anticoagulation Clinical Trial (PROACT)

## FDA Clinical Results

### Hemolysis

Postoperative serum lactate dehydrogenase (SLDH) levels are in the normal range, indicating minimal hemolysis<sup>8</sup>.

**Upper normal = 250**

<b>Valve</b>	<b>3-6 Months</b>	<b>1 Year</b>	<b>2 Years</b>
<b>Aortic</b>	222	225	229
<b>Mitral</b>	262	253	261
<b>Double</b>	238	246	243

## Morbid Event Rates

FDA clinical trials for the On-X valve produced the lowest morbid event rates of all carbon valves. The following tables compare On-X clinical data with those of other valves.

### On-X Morbid Event Rates Comparison – Aortic Position (% per pt-year)

Event	On-X <sup>1</sup>	SJM <sup>2</sup>	CMI <sup>3</sup>	ATS <sup>4</sup>	MH <sup>2</sup>
Thromboembolism	1.03	1.71	1.11	2.08	1.50
Thrombosis	0	0.19	0.03	0.00	0.20
Hemorrhage	0.59	2.00	1.82	1.96	0.80
Totals	1.62	3.90	2.96	4.04	2.50

### On-X Morbid Event Rates Comparison – Mitral Position (% per pt-year)

Event	On-X <sup>1</sup>	SJM <sup>2</sup>	CMI <sup>3</sup>	ATS <sup>4</sup>	MH <sup>2</sup>
Thromboembolism	1.00	1.86	2.05	4.00	2.00
Thrombosis	0	0.14	0.47	0.53	0.40
Hemorrhage	0.50	1.34	1.92	0.53	1.90
Totals	1.50	3.34	4.44	5.07	4.30

### Comparison On-X Morbid Event Rates – On-X Aortic vs. Tissue Aortic (% per pt-year)

Event	On-X <sup>1</sup> N=184 (2.8) 514**	SJM <sup>2</sup> N=577 (2.0) 1325**	CMI <sup>3</sup> N=868 (1.4) 1241**	ATS <sup>4</sup> N=267 (9.1) 2335**
Structural Dysfunction	0.0	0.0	0.0	2.0
Valve-Related Death	0.2	0.5	1.3	1.6
Thromboembolism	1.0	1.5	2.1	1.6
Totals	1.8	2.3	4.0	5.6

\* Combined subcoronary, full root and root inclusion.

\*\* Number of patients (average follow-up years) total patient years.

## References

1. Summary of Safety and Effectiveness, On-X<sup>®</sup> Prosthetic Heart Valve, Food and Drug Administration PMA P000037, May 30, 2001; P000037/S1 March 6, 2002; Primary European Trial P000037/R2, May 31, 2003.
2. Weighted averages calculated from references listed in the [On-X Experience, Volume 2, Issue 1, Sept. 2000. \(PDF\)](#)
3. Summary of Safety and Effectiveness, Carbomedics<sup>®</sup> Prosthetic Heart Valve, Food and Drug Administration PMA P900060, September 29, 1993.
4. Summary of Safety and Effectiveness, ATS Open Pivot<sup>®</sup> Bileaflet Heart Valve, Food and Drug Administration PMA P990046, October 13, 2000.
5. Summary of Safety and Effectiveness, St. Jude Medical Toronto SPV<sup>®</sup> Valve, Food and Drug Administration PMA P970030, November 26, 1997.
6. Summary of Safety and Effectiveness, Medtronic Freestyle<sup>®</sup> Aortic Root Bioprosthesis, Food and Drug Administration PMA P970031, November 24, 1997.
7. Frater et al., Long-term durability and patient functional status for the Carpentier-Edwards Perimount pericardial bioprosthesis in the aortic position, J Heart Valve Dis 1998;7:48-53.
8. Birnbaum D., Laczkovics A, Heidt M, Oelert H, Laufer G, Greve H, Pomar JL, Mohr F, Haverich A, Regensburger D., Examination of hemolytic potential with the On-X<sup>®</sup> prosthetic heart valve. J Heart Valve Dis. 2000; 9(1):142-5.